


One Shot Head Adaptor



- The One Shot Head Adaptor incorporates our unique disruption head and hydraulic operating and control system
- Quick and easy to use
- Allows loading of small fluid, viscous or non-flowing samples
- Contained during operation
- Available as an adaptor for all continuous flow models
- Collection cups provided can be transferred to and from ice or freezer for temperature control
- Easily dismantled for manual cleaning and autoclaving
- Materials in contact with sample, 316L stainless steel or equivalent, synthetic sapphire, EPDM, GLFPTFE & PEEF 450G

Maximum Pressure	40 KPSI (2700 bar)	30 KPSI (2000 bar)
Volume per Shot	1 - 10ml	1.5 - 13ml